

Easy

Safe

Fast



Animal Total RNA Isolation Kit
动物组织总RNA提取试剂盒



专注生命科学 · 让实验变得更简单



Foregene Co.,Ltd.

产品订购热线: 400-8566189

地址: 成都市温江区八一路北段18号三医创新中心A区D2栋

邮箱: info@foregene.com

网址: www.foregene.com

产品简介 INTRODUCTION

福际的动物组织RNA试剂盒采用本公司研制的离心柱和配方,可以从各种动物组织中高效率的提取得到高纯度高质量的总RNA。
提供高效的DNA-Cleaning Column,能轻松的让上清液和组织裂解物分离并吸附除去基因组DNA,操作简便、省时;
RNA-Only Column能高效的结合RNA,搭配独特的配方,可以同时处理大量样品。



无需低温冰浴

无需氯仿

无需DNase I


订购信息 ORDERING INFORMATION

| 产品名称 | 规格 | 货号 |
|---|------|----------|
| 动物组织总RNA提取试剂盒 Animal Total RNA Isolation Kit | 50T | RE-03011 |
| | 200T | RE-03014 |


下游应用 APPLICATION

动物组织总RNA提取试剂盒提取得到的总RNA可用于各种下游分子实验,例如:cDNA合成,RT-PCR、Real Time PCR、Northern Blot、Dot Blot、体外翻译、芯片分析、PolyA筛选、分子克隆和RNase保护分析等

实验流程 EXPERIMENTAL PROCESS



实验图例 EXPERIMENTAL DIAGRAM



Animal Total RNA Isolation Kit处理20mg小鼠新鲜样本,取5%纯化的总RNA 1%琼脂糖凝胶电泳。

引用文献 REFERENCES

- ◆ Zheng Q, Qin F, Luo R, et al. mRNA-Loaded Lipid-Like Nanoparticles for Liver Base Editing Via the Optimization of Central Composite Design. *Adv. Funct. Mater.* 2021, 31, 2011068. doi:10.1002/adfm.202011068. **IF:18.808(Foregene #.RE-03014)**
- ◆ He X, Hong W, Yang J, et al. Spontaneous apoptosis of cells in therapeutic stem cell preparation exert immunomodulatory effects through release of phosphatidylserine. *Signal Transduct Target Ther.* 2021 Jul 14;6(1):270. doi: 10.1038/s41392-021-00688-z. **IF:18.187(Foregene #.RE-03014)**
- ◆ Dai Z, Liu H, Liao J, et al. N7-Methylguanosine tRNA modification enhances oncogenic mRNA translation and promotes intrahepatic cholangiocarcinoma progression. *Mol Cell.* 2021 Jul 29;S1097-2765(21)00555-4. doi: 10.1016/j.molcel.2021.07.003. **IF:17.97(Foregene #.RE-03014)**
- ◆ Cao X, Shu Y, Chen Y, et al. Mettl14-Mediated m6A Modification Facilitates Liver Regeneration by Maintaining Endoplasmic Reticulum Homeostasis. *Cell Mol Gastroenterol Hepatol.* 2021;12(2):633-651. doi: 10.1016/j.jcmgh.2021.04.001. **IF:9.225(Foregene #.RE-03014)**

相关产品 RELATED PRODUCTS

第一链cDNA合成预混体系 Master Premix (For qPCR)

去基因组第一链cDNA合成试剂盒 (For qPCR) (With gDNase)

荧光定量PCR试剂盒-SYBR Green I 染料法

鼠尾直接PCR试剂盒—防PCR产物污染体系

用过都说好!

Make the experiment easier